

**Cleaning
And
Disinfection of
non- critical
Medical
Equipment**

Equipment / Site		Recommended way of how level disinfection and frequency
1	Arterial Blood gas machine	Use quaternary ammonia disinfectant wipes to clean wherever blood splashes occur and routinely 2x daily.
2	Breast Pump	<u>Non-Infected pt:</u> clean it c quaternary ammonia-based wipes <u>Infected pt.:</u> Dedicate one for each pt.'s mother upon discharge clean it c quaternary ammonia-based wipes.
3	Doppler probe	Clean it with quaternary ammonia based wipes after each use.

4	Blood pressure cuffs	<u>Non-Infected pt:</u> Clean it c alcohol wipes if c body fluid spillage sent to CSSD. <u>Infected pt:</u> Send to CSSD.
5	Chest tube clamp	Use alcohol-based wipes to clean after each use
6	Computers, Photocopy and Printers	Clean it c quaternary ammonia-based wipes daily
7	Crash cart trolley	Daily cleaning with detergent and water or ultrasol. If blood/body spillages to clean with 1000ppm sodium hypochlorite
8	Defibrillators	Clean it with quaternary ammonia-based wipes daily and after each use
9	Electrocardiography (ECG) machine	Prefer to use disposable ECG lead stickers. If not available to clean the attached metal leads with alcohol-based wipes between patients.
10	Blood Glucometer	Clean it with quaternary ammonia-based wipes after each use.

11	Humidifiers	<u>Non-infected patient:</u> Clean with detergent and water daily and dry it. The humidification water should be changed daily. <u>Infected patient:</u> send to CSSD
12	Incubators	If in use, clean with quaternary ammonia-based wipes once per shift. After discharge: <u>Non-infected patient:</u> add the step of cleaning with 1000 ppm sodium hypochlorite. <u>Infected patient:</u> Clean with detergent and water then clean with 1000 ppm sodium hypochlorite.
13	IV stand	Clean with detergent and water or ultrasol daily and dry it.
14	Laryngoscopes blade	Send to CSSD for sterilization
15	Laryngoscopes handle	Use alcohol-based wipes to clean the surface.
16	Medication trolley	Use alcohol-based wipes to clean the surface.
17	Monitors (vital signs, cardiac)	Clean it with quaternary ammonia alcohol-free based wipes.
18	Ophthalmoscopes	Use alcohol-based wipes to clean the surface.
19	Phototherapy capsule	If in use, clean with quaternary ammonia-based wipes once per shift. After discharge: <u>Non-infected patient:</u> add the step of cleaning with 1000 ppm sodium hypochlorite <u>Infected patient:</u> Clean with detergent and water then clean with 1000 ppm sodium hypochlorite.
20	Radiant warmers	Clean it with quaternary ammonia-based wipes after each use.
21	Scissors	Send to CSSD for sterilization
22	Stethoscopes	Clean with 70% alcohol swab after each patient.
23	Stretchers	<u>Non-infected patient:</u> Clean with detergent and water or ultrasol daily. <u>Infected patient:</u> Clean with detergent and water then clean with 1000ppm sodium hypochlorite.
24	Torch	Clean with 70% alcohol swab after each patient

ROYAL HOSPITAL
PROPER CLEANING AND WASHING OF
THE VIDEO LARYNGOSCOPE (C-MAC)

1. First, remove the black connection which is connected to the lens while using, connect the three-way adapter to the lens (to protect the lens).
2. Wash it with distilled water holding on to handle (ASH color) and let water run through.
3. Take 3-7.5ml of 3E-zyme and dilute it with 1L distilled water.
4. Soak it for 10 minutes.
5. Using the specific brush provided, scrub using the diluted 3E-zyme solution.
6. After cleaning with 3E-zyme, wash it again with distilled water.
7. Then soak it again in CIDEX for 15 minutes (take enough to dip, no need for dilution).
8. After that, take out from CIDEX and wash it THOROUGHLY in distilled water and DRY it.

NOTE: NEVER USE ALCOHOL WIPES OR ALCOHOL. USE ONLY SOAP WIPES FOR CLEANING THE CABLE.

For electronic equipment use non-alcohol ammonia based disinfectant wipes in both critical and general areas as well as in terminal cleaning. Use Ultrasol if wipes are not available or short of supply provided that J cloth squeezed well to prevent water leakage inside the equipment.

1) Terminal cleaning

- Ultrasol 10ml + 1L of water (2x)

2) Clostridium

Ultrasol 10ml + 1l H2O 1st cleaning

Haz tab

2.5 ----- 1.5 L H2o  for the 2nd cleaning

Or

4 ----- 2.5 L H2o